


It is quality rather than quantity that matters.

BANDAGES
SPONGES
COTTON WOOL



History and uses for weaved gauze

Gauze was traditionally woven in Palestine and the English word is said to derive from the place name for Gaza (Arabic: Ghazza), a center of weaving in the region.

Despite a prohibition on trade with non-Christians from religious authorities in medieval Europe, a fine type of silk known as gazzatum was imported from Gaza as early as the 13th century. Though members of religious orders in Europe were forbidden to wear it, the fabric won a place for itself and emerged into modern life as gauze.

According to the French government's online etymology

dictionary, the English form of the word derived from the French gaze, whose ultimate origin is uncertain, but is often attributed to the Arabic and Persian word gazz meaning "raw silk", which itself derived from the name of Gaza.[4] The same source says the existence of "an ancient textile industry in Gaza is not assured." and it is not known how the word entered into widespread use in European languages, with examples of first usages cited being the medieval Latin forms garza in Bologna in 1250 and gazzatum in Budapest in 1279.

Uses & types

Gauze was originally made of silk and was used for clothing. It is now used for many different things, including gauze sponges for medical purposes. When used as a medical dressing, gauze is generally made of cotton. It is especially useful for dressing wounds where other fabrics might stick to the burn or laceration. Many modern medical gauzes are covered with a plastic porous film such as Telfa or a polyblend which prevents direct contact and further minimizes wound adhesion. Also, it can be impregnated with a thick, creamy mixture of zinc oxide and calamine to promote healing, as in Unna's boot.

There are two types of wire gauze used in chemistry: a normally woven wire gauze, and one with a white circle imprinted on it. The latter has a ceramic composite intended to aid in the dispersal in heat. The normal set-up is a ring stand which a support ring attached to it. The wire gauze is centered on the ring stand as to allow the open flame to be in contact with it to a certain extent, depending on the desired heat. In the context of laboratory experimentation, the combination of the wire gauze plus the added ceramic composite greatly aids in the even dispersal of heat. This means that an object placed over this type of wire gauze will receive a more uniform heat from a Bunsen burner than the naked flame of it. This property is especially important in various chemical processes where precise, even heat dispersal is sought after, as a consistent heat is crucial for the successful maintenance of the chemical process intended. Because of this desired heat dispersal ability, wire gauze can also be used as a cooling surface for beakers, Erlenmeyer flasks or other vessels without damaging the supporting surface (workbench, tabletop, etc.)

PRODUCT DETAIL

1- GAUZES

2- BANDAGES

3- SPONGES

4- COTTON WOOL

BANDAGES

Bandages



SIZE & DISCRIPTION

```
ADHESIVE BANDAGE - 4.5 m x 6 cm
ADHESIVE BANDAGE - 4.5 m x 8 cm
COHESIVE BANDAGES - 4 m x 10 cm
COHESIVE BANDAGES - 4 m x 12 cm
COHESIVE BANDAGES - 4 m x 6 cm
COHESIVE BANDAGES - 4 m x 8 cm
ELASTIC COHESIVE BANDAGE - 20 m x 10 cm
ELASTIC COHESIVE BANDAGE - 20 m x 6 cm
ELASTIC COHESIVE BANDAGE - 20 m x 8 cm
ELASTIC COHESIVE BANDAGE - 4 m x 10 cm
ELASTIC COHESIVE BANDAGE - 4 m x 10 cm - blue
ELASTIC COHESIVE BANDAGE - 4 m x 10 cm - green
ELASTIC COHESIVE BANDAGE - 4 m x 10 cm - mixed
colours
ELASTIC COHESIVE BANDAGE - 4 m x 10 cm - red
ELASTIC COHESIVE BANDAGE - 4 m x 12 cm
ELASTIC COHESIVE BANDAGE - 4 m x 6 cm
ELASTIC COHESIVE BANDAGE - 4 m x 6 cm - blue
ELASTIC COHESIVE BANDAGE - 4 m x 6 cm - mixed
colours
ELASTIC COHESIVE BANDAGE - 4 m x 6 cm - red
ELASTIC COHESIVE BANDAGE - 4 m x 8 cm
ELASTIC COHESIVE BANDAGE - 4 m x 8 cm - blue
ELASTIC COHESIVE BANDAGE - 4 m x 8 cm - green
ELASTIC COHESIVE BANDAGE - 4 m x 8 cm - mixed
ELASTIC COHESIVE BANDAGE - 4 m x 8 cm - red
ELASTIC COHESIVE BANDAGE - 4 m x 8 cm -
yellow
ELASTICATED BANDAGE - 4.5 m x 10 cm
ELASTICATED BANDAGE - 4.5 m x 5 cm
```

GAUZES

Gauze Swabs



SIZE & DISCRIPTION

DESCRIPTION	Q.ty per box
Gauze Swabs - 11/7 - 40/40 - 5x5 - 8 plies without x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 7,5x7,5 - 8 plies without x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 7,5x10 - 8 plies without x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 7x9 - 8 plies without x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 10x10 - 8 plies without x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 10x20 - 8 plies without x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 7,5x7,5 - 12 plies without x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 10x10 - 12 plies without x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 10x10 - 12 plies with x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 10x20 - 12 plies without x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 5x5 - 16 plies without x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 5x10 - 16 plies without x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 7x9 - 16 plies without x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 7,5x7,5 - 16 plies without x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 10x10 - 16 plies without x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 10x20 - 16 plies without x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 5x5 - 8 plies with woven x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 7,5x7,5 - 8 plies with woven x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 7x9 - 8 plies with woven x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 10x10 - 8 plies with woven x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 10x20 - 8 plies with woven pvc x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 10x20 - 12 plies with x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 5x5 - 16 plies -with woven x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 7x9 - 16 plies with woven x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 7,5x7,5 - 16 plies with woven x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 10x20 - 16 plies with woven x-ray thread	100
Gauze Swabs - 11/7 - 40/40 - 10x10 - 16 plies with woven x-ray thread	100
Gauze swabs - 11/7 - 40/40 - 10x20 - 24 plies without x-ray thread	50
Gauze Swabs - 11/7 - 40/40 - 10x10 - 24 plies with woven x-ray thread	50
Gauze Swabs - 11/7 - 40/40 - 10x10 - 24 plies without x-ray thread	50

SPONGES

SPONGES





SIZE & DISCRIPTION

lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 10x60 cm - no loop	100
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 20x20 cm - no loop	100
lap.s ponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 25x25 cm - no loop	100
lap.s ponges 12/8 - 40/40 - with woven x-ray thread in PE- 4 plies - cm 30x30 cm no loop	100
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 30x45 cm no loop	100
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 40x40 cm no loop	100
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 45x45 cm - no loop	100
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 45x70 cm no loop	50
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 50x50 cm - no loop	100
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 10x60 cm - loop	100
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 20x20 cm - loop	100
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 25x25 cm - loop	100
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 30x30 cm - loop	100
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 30x45 cm - loop	100
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 40x40 cm - loop	100
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 45x45 cm - loop	100
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 45x70 cm - loop	50
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 45x70 cm - loop green color	50
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 50x50 cm - loop	100
lap.sponges 12/8 - 40/40 - with central woven x-ray thread in PE - 8 plies - cm 10x60 cm - no loop	100
lap.sponges 12/8 - 40/40 - with central woven x-ray thread in PE - 8 plies - cm 20x20 cm - no loop	100
lap.sponges 12/8 - 40/40 - with central woven x-ray thread in PE - 8 plies - cm 40x40 cm - no loop	100
lap.s ponges 12/12 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 20x20 cm - no loop	100
lap.sponges 12/12 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 25x25 cm - no loop	100
lap.sponges 12/12 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 30x30 cm - no loop	100
lap.sponges 12/12 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 40x40 cm - no loop	100
lap.sponges 12/12 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 50x50 cm - no loop	100
ap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 20x20 cm - loop 20 cm - Folded	1160/1 - tv
lap.sponges 12/12 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 20x20 cm - loop	100
lap.sponges 12/12 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 25x25 cm - loop	100
lap.sponges 12/12 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 30x30 cm - loop	100
lap.sponges 12/12 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 40x40 cm - loop	100
lap.sponges 12/12 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 50x50 cm - loop	100
lap.sponges 12/8 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 45x45 cm - loop	100
	100
lap.sponges 12/12 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 30x45 cm - loop	
lap.sponges 12/12 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 30x45 cm - loop lap.sponges 12/12 - 40/40 - with woven x-ray thread in PE - 4 plies - cm 30x45 cm - no loop	100
	100 50

COTTON WOOL

Cotton Wool



SIZE & DISCRIPTION

100% cotton interior, soft, pliable, non-linting, non-irritating, highlyabsorbent up to 6 times their weight. Will not adhere to mucous membranes, maintains shape better when wet, non sterile.

Carton of 10 boxes of 1000 pcs

- ✓ Cotton wool 250 g box of 20
- ✓ Cotton wool 1000 a box of 10
- ✓Cotton wool 100 g box of 50





18-Z Block, 3rd Floor Commercial Area, III, DHA, Lahore-Pakistan

Tel: 0092 42 35693147 Fax: 0092 42 35693137 Cell: 0092 0303 444 2020

Info@ccintgroup.com s.k@ccintgroup.com

